

For Research Use Only

UFL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83047-4-PBS



Basic Information

Catalog Number:	83047-4-PBS	GenBank Accession Number:	BC036379	Purification Method:	Protein A purification
Source:	Rabbit	GeneID (NCBI):	23376	CloneNo.:	230363F10
Isotype:	IgG	UNIPROT ID:	O94874		
Immunogen Catalog Number:	AG23362	Full Name:	KIAA0776		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00135-1, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-3-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00135-2, KIAA0776 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83047-1-PBS. Detection antibody: 83047-4-PBS. Standard: Ag23362. Range: 0.313-20 ng/mL.