

For Research Use Only

KIAA0776 Recombinant antibody

Catalog Number: 83047-6-RR



Basic Information

Catalog Number:

83047-6-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23362

GenBank Accession Number:

BC036379

GeneID (NCBI):

23376

UNIPROT ID:

O94874

Full Name:

KIAA0776

Purification Method:

Protein A purification

CloneNo.:

230363A8

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

