

For Research Use Only

# TM9SF4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83048-1-PBS

## Basic Information

Catalog Number:	83048-1-PBS	GenBank Accession Number:	BC021107	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	9777	CloneNo.:	230344C12
Source:	Rabbit	UNIPROT ID:	Q92544		
Isotype:	IgG	Full Name:	transmembrane 9 superfamily protein member 4		
Immunogen Catalog Number:	AG31253	Calculated MW:	642 aa, 75 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00139-1, TM9SF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83048-2-PBS. Detection antibody: 83048-1-PBS. Standard: Ag31253. Range: 1.25-20 ng/mL