

For Research Use Only

TM9SF4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83048-2-PBS

Basic Information

Catalog Number:	83048-2-PBS	GenBank Accession Number:	BC021107	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	9777	CloneNo.:	230344D8
Source:	Rabbit	UNIPROT ID:	Q92544		
Isotype:	IgG	Full Name:	transmembrane 9 superfamily protein member 4		
Immunogen Catalog Number:	AG31253	Calculated MW:	642 aa, 75 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00139-1, TM9SF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83048-2-PBS. Detection antibody: 83048-1-PBS. Standard: Ag31253. Range: 1.25-20 ng/mL