

For Research Use Only



# Integrin beta 3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83053-1-PBS

## Basic Information

Catalog Number:	83053-1-PBS	GenBank Accession Number:	BC127667	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	3690	CloneNo.:	230292B2
Source:	Rabbit	UNIPROT ID:	P05106		
Isotype:	IgG	Full Name:	integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)		
		Calculated MW:	87 kDa		

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00047-2, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-1-PBS. Standard: SY00181. Range: 0.312-5 ng/mL