

For Research Use Only

Integrin beta 3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83053-2-PBS

Basic Information

Catalog Number:

83053-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC127667

GeneID (NCBI):

3690

UNIPROT ID:

P05106

Full Name:

integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)

Calculated MW:

87 kDa

Purification Method:

Protein A purification

CloneNo.:

230292B6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00047-1, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-3-PBS. Standard: SY00181. Range: 0.312-5 ng/mL



Cytometric bead array standard curve ofMP00047-2, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-1-PBS. Standard: SY00181. Range: 0.312-5 ng/mL