

For Research Use Only



Integrin beta 3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83053-3-PBS

Basic Information

Catalog Number:	83053-3-PBS	GenBank Accession Number:	BC127667	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	3690	CloneNo.:	230292F12
Source:	Rabbit	UNIPROT ID:	P05106		
Isotype:	IgG	Full Name:	integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)		
		Calculated MW:	87 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00047-1, Integrin beta-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83053-2-PBS. Detection antibody: 83053-3-PBS. Standard: SY00181. Range: 0.312-5 ng/mL.