

For Research Use Only

GABRA2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83057-1-PBS



Basic Information

Catalog Number:

83057-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29301

GenBank Accession Number:

NM_000807

GeneID (NCBI):

2555

UNIPROT ID:

P47869

Full Name:

gamma-aminobutyric acid (GABA) A receptor, alpha 2

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

230158C12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00067-1, GABRA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83057-3-PBS. Detection antibody: 83057-1-PBS. Standard: Ag29301. Range: 0.312-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83057-1-PBS against Human GABRA2 was performed. The affinity constant is 64.0 pM.