

For Research Use Only

SOX2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83059-3-PBS



Basic Information

Catalog Number:	83059-3-PBS	GenBank Accession Number:	BC013923	Purification Method:	Protein A purification
Concentration:	1000 µg/ml	GeneID (NCBI):	6657	CloneNo.:	230225E7
Source:	Rabbit	UNIPROT ID:	P48431		
Isotype:	IgG	Full Name:	SRV (sex determining region Y)-box 2		
Immunogen Catalog Number:	AG0612	Calculated MW:	34 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00077-1, SOX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83059-1-PBS. Detection antibody: 83059-3-PBS. Standard: Ag0612. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00077-2, SOX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83059-2-PBS. Detection antibody: 83059-3-PBS. Standard: Ag0612. Range: 0.625-20 ng/mL