

For Research Use Only

ABCA13 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83060-2-PBS



Basic Information

Catalog Number:

83060-2-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_152701

GeneID (NCBI):

154664

UNIPROT ID:

Q86UQ4

Full Name:

ATP-binding cassette, sub-family A (ABC1), member 13

Calculated MW:

576 kDa

Purification Method:

Protein A purification

CloneNo.:

230288G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00084-1, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-1-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL.



Cytometric bead array standard curve of MP00084-2, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-3-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL.