

For Research Use Only

# ABCA13 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83060-3-PBS

## Basic Information

Catalog Number:

83060-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_152701

GeneID (NCBI):

154664

UNIPROT ID:

Q86UQ4

Full Name:

ATP-binding cassette, sub-family A (ABC1), member 13

Calculated MW:

576 kDa

Purification Method:

Protein A purification

CloneNo.:

230288H4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00084-2, ABCA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83060-3-PBS. Detection antibody: 83060-2-PBS. Standard: SY00509. Range: 0.312-20 ng/mL.