

For Research Use Only

TRPS1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83062-2-PBS

Basic Information

Catalog Number: 83062-2-PBS	GenBank Accession Number: BC125020	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7227	CloneNo.: 230051F12
Source: Rabbit	UNIPROT ID: Q9UHF7	Affinity: $K_D = 4.73 \times 10^{-11}$
Isotype: IgG	Full Name: trichorhinophalangeal syndrome I	$K_{off} = 1.63 \times 10^{-5}$
Immunogen Catalog Number: AG33959	Calculated MW: 1294 aa, 143 kDa	$K_{on} = 3.44 \times 10^5$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00096-2, TRPS1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83062-2-PBS. Detection antibody: 83062-1-PBS. Standard: Ag33959. Range: 1.25-20 ng/mL.