

For Research Use Only

MOG Recombinant antibody, PBS Only (Capture)

Catalog Number: 83063-1-PBS



Basic Information

Catalog Number:

83063-1-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3273

GenBank Accession Number:

BC035938

GeneID (NCBI):

4340

UNIPROT ID:

Q16653

Full Name:

myelin oligodendrocyte glycoprotein

Calculated MW:

295 aa, 34 kDa

Purification Method:

Protein A purification

CloneNo.:

230291C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00110-1, MOG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83063-1-PBS. Detection antibody: 83063-4-PBS. Standard: Ag3273. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00110-3, MOG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83063-1-PBS. Detection antibody: 83063-3-PBS. Standard: Ag3273. Range: 0.625-20 ng/mL