

For Research Use Only

# Mouse Tnfrsf11b Recombinant antibody, PBS Only (Detector)



Catalog Number: 83065-4-PBS

## Basic Information

Catalog Number:

83065-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC049782.1

GeneID (NCBI):

18383

UNIPROT ID:

Q6PI12

Full Name:

tumor necrosis factor receptor  
superfamily, member 11b  
(osteoprotegerin)

Calculated MW:

46KD

Purification Method:

Protein A purification

CloneNo.:

230325G4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00123-1, Mouse Tnfrsf11b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83065-1-PBS. Detection antibody: 83065-4-PBS. Standard: Eg0598. Range: 0.625-20 ng/mL.