

For Research Use Only

RNF13 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83066-2-PBS

Basic Information

Catalog Number: 83066-2-PBS	GenBank Accession Number: BC009803	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 11342	CloneNo.: 230351D11
Source: Rabbit	UNIPROT ID: O43567	
Isotype: IgG	Full Name: ring finger protein 13	
Immunogen Catalog Number: AG8987	Calculated MW: 381 aa, 43 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00008-1, RNF13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83066-3-PBS. Detection antibody: 83066-2-PBS. Standard: Ag8987. Range: 0.625-10 ng/mL