

For Research Use Only

# ECM1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83067-3-PBS

## Basic Information

Catalog Number:

83067-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC023505

GeneID (NCBI):

1893

UNIPROT ID:

Q16610

Full Name:

extracellular matrix protein 1

Calculated MW:

540 aa, 61 kDa

Purification Method:

Protein A purification

CloneNo.:

230321G12

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00011-2, ECM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83067-3-PBS. Detection antibody: 83067-5-PBS. Standard: Eg0533. Range: 0.625-20 ng/mL