

For Research Use Only

IFNAR2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83068-1-PBS

Basic Information

Catalog Number:	83068-1-PBS	GenBank Accession Number:	BC002793	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	3455	CloneNo.:	23022485
Source:	Rabbit	UNIPROT ID:	P48551		
Isotype:	IgG	Full Name:	interferon (alpha, beta and omega) receptor 2		
Immunogen Catalog Number:	AG0191	Calculated MW:	331 aa, 37 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00016-2, IFNAR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83068-1-PBS. Detection antibody: 83068-2-PBS. Standard: Ag0191. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00016-3, IFNAR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83068-1-PBS. Detection antibody: 83068-2-PBS. Standard: Ag0191. Range: 0.625-20 ng/mL