

For Research Use Only

# IFNAR2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83068-3-PBS

## Basic Information

Catalog Number: 83068-3-PBS	GenBank Accession Number: BC002793	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3455	CloneNo.: 230224D11
Source: Rabbit	UNIPROT ID: P48551	
Isotype: IgG	Full Name: interferon (alpha, beta and omega) receptor 2	
Immunogen Catalog Number: AG0191	Calculated MW: 331 aa, 37 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve ofMP00016-1, IFNAR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83068-3-PBS. Detection antibody: 83068-4-PBS. Standard: Ag0191. Range: 0.625-20 ng/mL