

For Research Use Only

# IFNAR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83068-7-PBS



## Basic Information

Catalog Number: 83068-7-PBS	GenBank Accession Number: NM_207585.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3455	CloneNo.: 250313B10
Source: Rabbit	UNIPROT ID: P48551-1	
Isotype: IgG	Full Name: interferon (alpha, beta and omega) receptor 2	
Immunogen Catalog Number: EG2928	Calculated MW: 58 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean IFNAR2 concentration was determined to be 24.5 ng/mL in human serum.



Cytometric bead array standard curve of MP00016-4, IFNAR2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83068-7-PBS. Detection antibody: 83068-6-PBS. Standard: Eg2928. Range: 1.563-100 ng/mL.