

For Research Use Only

# ICAM-1/CD54 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83069-2-PBS



## Basic Information

Catalog Number:

83069-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC015969

GeneID (NCBI):

3383

ENSEMBL Gene ID:

ENSG00000090339

UNIPROT ID:

P05362

Full Name:

intercellular adhesion molecule 1

Calculated MW:

90 kDa

Purification Method:

Protein A purification

CloneNo.:

230342F2

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00028-4, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-5-PBS. Detection antibody: 83069-2-PBS. Standard: Eg0570. Range: 78.1-5000 pg/mL.



Sandwich ELISA standard curve of MP00028-5, ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83069-1-PBS. Detection antibody: 83069-2-PBS. Standard: Eg0570. Range: 39-2500 pg/mL.