

For Research Use Only

PARP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83071-2-PBS



Basic Information

Catalog Number:

83071-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005484

GeneID (NCBI):

10038

UNIPROT ID:

Q9UGN5

Full Name:

poly (ADP-ribose) polymerase 2

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

230386E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00045-3, PARP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83071-2-PBS. Detection antibody: 83071-3-PBS. Standard: SY01205. Range: 0.078-5 ng/mL.