

For Research Use Only

# PARP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83071-4-PBS



## Basic Information

Catalog Number:

83071-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005484

GeneID (NCBI):

10038

UNIPROT ID:

Q9UGN5

Full Name:

poly (ADP-ribose) polymerase 2

Calculated MW:

66 kDa

Purification Method:

Protein A purification

CloneNo.:

230386H5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00045-2, PARP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83071-1-PBS. Detection antibody: 83071-4-PBS. Standard: SY01205. Range: 0.078-5 ng/mL.