

For Research Use Only

MFRP Recombinant antibody, PBS Only (Detector)



Catalog Number: 83075-3-PBS

Basic Information

Catalog Number: 83075-3-PBS	GenBank Accession Number: BC160191	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 83552	CloneNo.: 230418D5
Source: Rabbit	UNIPROT ID: Q9BY79	Affinity: $K_D = 1.16 \times 10^{-9}$
Isotype: IgG	Full Name: membrane frizzled-related protein	$K_{off} = 3.62 \times 10^{-5}$ $K_{on} = 3.12 \times 10^4$
Immunogen Catalog Number: AG34302	Calculated MW: 62 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00066-1, MFRP Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83075-1-PBS. Detection antibody: 83075-3-PBS. Standard: Ag34302. Range: 0.625-10 ng/mL