

For Research Use Only

SLC4A9 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83076-1-PBS

Basic Information

Catalog Number:	83076-1-PBS	GenBank Accession Number:	NM_001258428	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	83697	CloneNo.:	230419A7
Source:	Rabbit	UNIPROT ID:	Q96Q91		
Isotype:	IgG	Full Name:	solute carrier family 4, sodium bicarbonate cotransporter, member 9		
Immunogen Catalog Number:	AG34492	Calculated MW:	983 aa, 108 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00073-2, SLC4A9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83076-3-PBS. Detection antibody: 83076-1-PBS. Standard: Ag34492. Range: 2.5-20 ng/mL.