

For Research Use Only

STT3A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83077-2-PBS

Basic Information

Catalog Number:

83077-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31994

GenBank Accession Number:

BC020965

GeneID (NCBI):

3703

UNIPROT ID:

P46977

Full Name:

STT3, subunit of the oligosaccharyltransferase complex, homolog A (*S. cerevisiae*)

Calculated MW:

705 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

230426D2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00080-1, STT3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83077-3-PBS. Detection antibody: 83077-2-PBS. Standard: Ag31994. Range: 0.625-20 ng/mL