

For Research Use Only

p300 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83078-3-PBS



Basic Information

Catalog Number:

83078-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001429

GeneID (NCBI):

2033

UNIPROT ID:

Q09472

Full Name:

E1A binding protein p300

Calculated MW:

264 kDa

Purification Method:

Protein A purification

CloneNo.:

230369H1

Affinity:

$K_D = 1.22 \times 10^{-9}$

$K_{off} = 1.15 \times 10^{-4}$

$K_{on} = 9.45 \times 10^4$

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00310-2, p300 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83078-4-PBS. Detection antibody: 83078-3-PBS. Standard: SY00911. Range: 0.313-5 ng/mL.