

For Research Use Only

CA11 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83082-3-PBS

Basic Information

Catalog Number: 83082-3-PBS	GenBank Accession Number: BC002662	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 770	CloneNo.: 230390E5
Source: Rabbit	UNIPROT ID: O75493	
Isotype: IgG	Full Name: carbonic anhydrase XI	
Immunogen Catalog Number: AG33373	Calculated MW: 36 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00109-1, CA11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83082-3-PBS. Detection antibody: 83082-2-PBS. Standard: Ag33373. Range: 1.25-20 ng/mL



Cytometric bead array standard curve ofMP00109-2, CA11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83082-1-PBS. Detection antibody: 83082-3-PBS. Standard: Ag33373. Range: 1.25-20 ng/mL