

For Research Use Only

ABCA6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83084-1-PBS

Basic Information

Catalog Number: 83084-1-PBS	GenBank Accession Number: NM_080284	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 23460	CloneNo.: 230223A11
Source: Rabbit	UNIPROT ID: Q8N139	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family A (ABC1), member 6	
Immunogen Catalog Number: AG29524	Calculated MW: 184 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00118-1, ABCA6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83084-1-PBS. Detection antibody: 83084-2-PBS. Standard: Ag29524. Range: 2.5-20 ng/mL.