

For Research Use Only

ABCA6 Recombinant antibody, PBS Only



Catalog Number: 83084-4-PBS

Basic Information

Catalog Number: 83084-4-PBS	GenBank Accession Number: NM_080284	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 23460	CloneNo.: 230223E8
Source: Rabbit	UNIPROT ID: Q8N139	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family A (ABC1), member 6	
Immunogen Catalog Number: AG29524	Calculated MW: 184 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
Human

Background Information

ABCA6, a member of the superfamily of ATP-binding cassette (ABC) transporters, is involved in the energy-dependent transport of a wide variety of substrates across membranes. The human ABCA6 is ubiquitously expressed in the organism, with high levels in liver, lung, heart, and brain(PMID: 11478798).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83084-4-RR against Human ABCA6 were performed. The affinity constant is 11.6 nM.