

For Research Use Only

# SLC35F4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83085-2-PBS

## Basic Information

Catalog Number:

83085-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34518

GenBank Accession Number:

NM\_001206920

GeneID (NCBI):

341880

Full Name:

solute carrier family 35, member F4

Calculated MW:

521 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

230360H1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00122-1, SLC35F4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83085-1-PBS. Detection antibody: 83085-2-PBS. Standard: Ag34518. Range: 1.25-20 ng/mL