

For Research Use Only

MOGAT2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83086-2-PBS



Basic Information

Catalog Number:

83086-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15806

GenBank Accession Number:

BC103878

GeneID (NCBI):

80168

UNIPROT ID:

Q3SYC2

Full Name:

monoacylglycerol O-acyltransferase

2

Calculated MW:

334 aa, 38 kDa

Purification Method:

Protein A purification

CloneNo.:

230365F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00124-1, MOGAT2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83086-1-PBS. Detection antibody: 83086-2-PBS. Standard: Ag15806. Range: 0.625-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83086-2-PBS against Human MOGAT2 was performed. The affinity constant is 1.13 nM.