

For Research Use Only

GDNF Recombinant antibody, PBS Only



Catalog Number: 83087-2-PBS

Basic Information

Catalog Number:

83087-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24068

GenBank Accession Number:

BC069369

GeneID (NCBI):

2668

UNIPROT ID:

P39905

Full Name:

glial cell derived neurotrophic factor

Calculated MW:

211 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

230358G6

Applications

Tested Applications:

ELISA, WB

Species Specificity:

mouse, human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 83087-2-RR (GDNF antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83087-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83087-2-RR against Human GDNF were performed. The affinity constant is below 1 pM.