

For Research Use Only

PRSS16 Recombinant antibody

Catalog Number: 83091-1-RR



Basic Information

Catalog Number:

83091-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33961

GenBank Accession Number:

BC156928

GeneID (NCBI):

10279

UNIPROT ID:

Q9NQE7

Full Name:

protease, serine, 16 (thymus)

Purification Method:

Protein A purification

CloneNo.:

230190D11

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Background Information

PRSS16, also referred to as Thymic Specific Serine Protease (TSSP) is specifically expressed by cortical epithelial cells of the thymus and localizes around the thymic capsule and vessels (PMID: 12144812).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83091-1-RR against Human PRSS16 were performed. The affinity constant is below 1 pM.