

For Research Use Only

# PAPD5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83098-4-PBS

## Basic Information

Catalog Number:

83098-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001040285

GeneID (NCBI):

64282

UNIPROT ID:

Q8NDF8

Full Name:

PAP associated domain containing 5

Calculated MW:

64 kDa

Purification Method:

Protein A purification

CloneNo.:

230368D8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00145-1, PAPD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83098-2-PBS. Detection antibody: 83098-4-PBS. Standard: SY01262. Range: 0.313-10 ng/mL



Cytometric bead array standard curve ofMP00145-2, PAPD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83098-3-PBS. Detection antibody: 83098-4-PBS. Standard: SY01262. Range: 0.313-20 ng/mL