

For Research Use Only

TGOLN2,TGN46 Recombinant antibody, PBS Only (Capture)

Catalog Number:83099-4-PBS



Basic Information

Catalog Number:

83099-4-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4470

GenBank Accession Number:

BC028219

GeneID (NCBI):

10618

UNIPROT ID:

O43493

Full Name:

trans-golgi network protein 2

Calculated MW:

447 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

230376D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00152-3, TGN46 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83099-4-PBS. Detection antibody: 83099-3-PBS. Standard: Ag4470. Range: 0.625-20 ng/mL