

For Research Use Only

TNFR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83101-3-PBS



Basic Information

Catalog Number:

83101-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230328E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00777-2, TNFR2/CD120b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83101-3-PBS. Detection antibody: 98010-1-PBS. Standard: Eg0086. Range: 0.625-20 ng/mL.