

For Research Use Only

# TNFR2/CD120b Recombinant antibody, PBS Only

Catalog Number: 83101-4-PBS



## Basic Information

Catalog Number:

83101-4-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0086

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor  
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230328F5

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83101-4-PBS against Human TNFR2 were performed. The affinity constant is 3.77 nM.