

For Research Use Only

# TNFR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83101-8-PBS



## Basic Information

Catalog Number:

83101-8-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0086

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor  
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230328B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

