

For Research Use Only

Mouse Lipocalin-2/NGAL Recombinant antibody, PBS Only (Capture)

Catalog Number: 83102-6-PBS



Basic Information

Catalog Number:

83102-6-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008491.1

GeneID (NCBI):

16819

UNIPROT ID:

P11672

Full Name:

lipocalin 2

Calculated MW:

23KD

Purification Method:

Protein A purification

CloneNo.:

230372D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00302-3, Mouse NGAL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83102-6-PBS. Detection antibody: 83102-4-PBS. Standard: Eg0572. Range: 0.625-20 ng/mL.