

For Research Use Only

DDX19A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83106-3-PBS

Basic Information

Catalog Number:

83106-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7362

GenBank Accession Number:

BC005162

GeneID (NCBI):

55308

UNIPROT ID:

Q9NUU7

Full Name:

DEAD (Asp-Glu-Ala-As) box polypeptide 19A

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

230232D3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00140-1, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-3-PBS. Detection antibody: 83106-6-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00140-2, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-3-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL