

For Research Use Only

DDX19A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83106-4-PBS

Basic Information

Catalog Number:	83106-4-PBS	GenBank Accession Number:	BC005162	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	55308	CloneNo.:	230232D7
Source:	Rabbit	UNIPROT ID:	Q9NUU7		
Isotype:	IgG	Full Name:	DEAD (Asp-Glu-Ala-As) box polypeptide 19A		
Immunogen Catalog Number:	AG7362	Calculated MW:	54 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00140-2, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-3-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00140-3, DDX19A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83106-5-PBS. Detection antibody: 83106-4-PBS. Standard: Ag7362. Range: 0.625-20 ng/mL