

For Research Use Only

# TJP2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 83111-3-PBS



## Basic Information

Catalog Number:

83111-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13534

GenBank Accession Number:

BC027592

GeneID (NCBI):

9414

UNIPROT ID:

Q9UDY2

Full Name:

tight junction protein 2 (zona occludens 2)

Calculated MW:

134 kDa

Purification Method:

Protein A purification

CloneNo.:

230474F8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00241-1, ZO-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83111-2-PBS. Detection antibody: 83111-3-PBS. Standard: Ag13534. Range: 0.625-20 ng/mL