

For Research Use Only

# CELSR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83115-3-PBS



## Basic Information

Catalog Number:

83115-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001408

GeneID (NCBI):

1952

UNIPROT ID:

Q9HCU4

Full Name:

cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)

Calculated MW:

318 kDa

Purification Method:

Protein A purification

CloneNo.:

230319D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00165-1, Human CELSR2  
Recombinant Matched Antibody Pair, PBS Only  
measured by Cytometric bead array. Capture  
antibody: 83115-4-PBS. Detection antibody: 83115-  
3-PBS. Standard:null. Range: null.