

For Research Use Only

# WHSC1L1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83116-4-PBS



## Basic Information

Catalog Number:

83116-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1894

GenBank Accession Number:

BC012059

GeneID (NCBI):

54904

UNIPROT ID:

Q9BZ95

Full Name:

Wolf-Hirschhorn syndrome candidate

1-like 1

Calculated MW:

1437 aa, 162 kDa

Purification Method:

Protein A purification

CloneNo.:

230499F12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00189-1, WHSC1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83116-6-PBS. Detection antibody: 83116-4-PBS. Standard: Ag1894. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00189-2, WHSC1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83116-2-PBS. Detection antibody: 83116-4-PBS. Standard: Ag1894. Range: 0.625-20 ng/mL