

For Research Use Only

# INO80 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83117-2-PBS



## Basic Information

Catalog Number:

83117-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC146785

GeneID (NCBI):

54617

UNIPROT ID:

Q9ULG1

Full Name:

INO80 homolog (S. cerevisiae)

Calculated MW:

1556 aa, 177 kDa

Purification Method:

Protein A purification

CloneNo.:

230440F5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00159-1, INO80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83117-3-PBS. Detection antibody: 83117-2-PBS. Standard: SY00250. Range: 0.078-2.5 ng/mL