

For Research Use Only

Frizzled 5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83119-2-PBS



Basic Information

Catalog Number:

83119-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16082

GenBank Accession Number:

BC114346

GeneID (NCBI):

7855

UNIPROT ID:

Q13467

Full Name:

frizzled homolog 5 (Drosophila)

Calculated MW:

585 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

230442B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00265-1, FZD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83119-4-PBS. Detection antibody: 83119-2-PBS. Standard: Ag16082. Range: 1.25-20 ng/mL.



Cytometric bead array standard curve of MP00265-2, FZD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83119-3-PBS. Detection antibody: 83119-2-PBS. Standard: Ag16082. Range: 1.25-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83119-2-PBS against Human Frizzled 5 was performed. The affinity constant is 8.29 nM.