

For Research Use Only

OPRM1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83120-4-PBS



Basic Information

Catalog Number:

83120-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24806

GenBank Accession Number:

BC074927

GeneID (NCBI):

4988

UNIPROT ID:

P35372

Full Name:

opioid receptor, mu 1

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

230444G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00271-1, OPRM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83120-4-PBS. Detection antibody: 83120-2-PBS. Standard: Ag24806. Range: 1.25-20 ng/mL.



Cytometric bead array standard curve of MP00271-2, OPRM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83120-3-PBS. Detection antibody: 83120-4-PBS. Standard: Ag24806. Range: 0.625-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 83120-4-PBS against Human OPRM1 was performed. The affinity constant is 11.9 pM.