

For Research Use Only

SLC25A28 Recombinant antibody, PBS Only (Capture)

Catalog Number:83127-1-PBS



Basic Information

Catalog Number:

83127-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34507

GenBank Accession Number:

BC076399

GeneID (NCBI):

81894

UNIPROT ID:

Q96A46

Full Name:

solute carrier family 25, member 28

Purification Method:

Protein A purification

CloneNo.:

230392A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00142-3, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-1-PBS. Detection antibody: 83127-3-PBS. Standard: Ag34507. Range: 0.625-20 ng/mL.