

For Research Use Only

SLC25A28 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83127-2-PBS

Basic Information

Catalog Number: 83127-2-PBS	GenBank Accession Number: BC076399	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 81894	CloneNo.: 230392C1
Source: Rabbit	UNIPROT ID: Q96A46	
Isotype: IgG	Full Name: solute carrier family 25, member 28	
Immunogen Catalog Number: AG34507		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00142-1, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-4-PBS. Detection antibody: 83127-2-PBS. Standard: Ag34507. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00142-2, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-2-PBS. Detection antibody: 83127-5-PBS. Standard: Ag34507. Range: 1.25-20 ng/mL