

For Research Use Only

# SLC25A28 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83127-4-PBS

## Basic Information

Catalog Number:

83127-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34507

GenBank Accession Number:

BC076399

GeneID (NCBI):

81894

UNIPROT ID:

Q96A46

Full Name:

solute carrier family 25, member 28

Purification Method:

Protein A purification

CloneNo.:

230392D4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00142-1, SLC25A28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83127-4-PBS. Detection antibody: 83127-2-PBS. Standard: Ag34507. Range: 0.625-20 ng/mL.